

Exercise at the Spectrometer

General

The script is written for Bruker spectrometers. The color code means, red for all commands and green for parameters in the TopSpin commando line. The script gives additional information to the recommended Bruker manuals: Beginner Guide, 1D and 2D Step-by-Step -Basic & - Advanced which can be found under the help button in TopSpin.

set up a 1D NMR experiment, general procedure:

- Temperature adjustment
- Load sample into the magnet
- create a new data set (e.g. *rpar PROTON all*)
- tuning the probe
- Iock & shimming
- determine the acquisition parameters
- start the experiment
- processing
- plotting

Probe temperature

Before you load a new sample into the magnet check the actual temperature in *edte*. For cryogenic probe heads it is necessary that the gas flow is 670 l/h and the temperature calibration is active.

Load the sample into the magnet

Use the Bruker tool at the spectrometer to center the sample in the spinner around the active volume. For cryogenic probe heads the maximum sample depth is 21mm. The right orientation is shown in Figure 1, step 5. In addition to the standard NMR tubes for Bio-NMR in a 5mm probe head the following tubes can be used:

Tubes	Diameter	Sample volume	Coasts	Comments
High quality	5mm	550 μl	10 – 20 €	
High quality	3mm	180 μl	10 – 20 €	for lossy biological sample
Shigemi tube	5mm	280 μl	110€	for expensive sample
Shape tube	oval tube	340 μl	260 €	for lossy biological sample





Figure1 from the Bruker manual beginner guide.

For Shigemi tubes the active volume between the bottom glass and the plunger glass should be centered. It is important to expel air bubbles from the sample, especially when Shigemi tubes are used.

Create a new data set

New or edc generate a new data set name. A dataset is defined by 5 parameters:

<dir>:</dir>	The directory where the dataset is located
<user>:</user>	The name of the TopSpin user who did the experiment
<name>:</name>	Identifier of the sample
<expno>:</expno>	Contains the raw data of one experiment (= FID)
<procno>:</procno>	Contains the processed data of one experiment (= spectrum)

The command rpar reads a parameter set (experiment) to the current dataset. When it is entered without arguments, rpar opens a dialog box with a list of the available parameter sets. The naming of the parameter sets corresponds to the name of the pulse sequences, where the first characters specify the type of experiment, e.g. DEPT, COSY, NOESY etc. Further properties of the pulseprogram are indicated by a two-character code, which is added to the name in alphabetical order. The Bruker naming convention can be checked in edcpul Pulprog.info. When using the parameter set to start NMR-experiments it is important to keep in mind that the pulseprogram is only one of the parameters.

Commands to change the data structure

- edc, new
 create a new data set
 - *search* find an existing data set
- *re expno* with this command you change the expnos
- *rep procno* with this command you change the procnos
- *wrp procno* copy the processed data to a new procno
- wrpa expno
 copy the FID and processed data to a new expno



tuning the probe

The probehead is a LC resonance circuit (like an old radio receiver) which needs to be adjusted to the frequency used. This procedure is called tuning & matching. While tuning means optimizing to the frequency, matching will minimize the reflected power.

Procedure for tuning and matching the probe head

define the channels in edasp



Set default paths

- wobb start the wobble process
- Color code of T + M screws:

0	yellow:	¹ H Proton
0	blue:	¹³ C Carbon
0	red:	¹⁵ N Nitrogen
0	grey:	² H Deuterium (only CryoProbe TM)

 change the tuning und matching screws at the probe head to reach a minimum at the red central line.





- Automated tuning and matching:
 - Special probe head required: ATM
 - Commands: atm (fully automated) or atmm (manually)
- Procedure T + M of ²H Deuterium (only for cryogenic Probes)
 - Generate a new data set
 - read parameter set for deuterium gradient shimming: rpar gradshim1d2h & wobb tune and match the probe
 - Afterwards change back to ¹H-1D-dataset and do ii for lock again (ii initializes the spectrometer. It is only a software command and can be used any time)

Lock and shimming

The magnetic field B₀ requires stabilization of the main field (LOCK) and control of the field homogeneity (SHIM).

LOCK

- continuously determines the frequency of ²H signal of the solvent (deuterated solvents)
- adds a small extra field to the main field of the magnet to keep the overall field constant
- ²H signal also used for shimming

SHIM

- additional coils outside the sample used for adjustment of the homogeneity of the B₀ field (e.g. Z, Z2, Z3, X, Y....)
- An NMR Spectrometer has <u>cryo shims</u> and <u>room temperature shims</u>

The lock channel can be thought as a *,complete independent second spectrometer within the main spectrometer'* (Larmor frequency: $\omega_0 = \gamma B_0$). An additional field H_0 is created: $\omega_0 = \gamma B_0$



 $\gamma(B_0+H_0)$. B_0 is not constant, but (B_0+H_0) can be held constant by control technology (field regulation). The lock signal is recorded by two quadrature channels in absorption and dispersion. The dispersion signal is used for field stabilisation.

LOCK Parameter

- LOCK gain The gain setting for the display only.
- LOCK Power

The ²H transmitter output power. Due to the different relaxation behavior of ²H for individual solvents, the lock power has to be adjusted for each solvent. Too high lock power will result in an unstable lock signal due the saturation. To check if the look power is close to saturation, change the lock power of -1dB and decrease lock gain by -1dB. If the lock line is at same position as before, repeat the procedure.

LOCK phase

If the lock phase is not adjusted correctly, absorption and dispersion signals will be mixed. Non-pure phases will result in imperfect field homogenisation (shimming) and field shifts during experiment using pulsed field gradients.

- Regulation parameters
 - Loop Gain: how *strong* to react on field disturbance
 - Loop Time: how *fast* to react on field disturbance
 - Loop Filter: *smoothing* the lock signal to remove noise, low pass filter

Wrong settings will result in unstable signal position of the lock, which can give suppression artifacts or additional noise around the signal in NMR-spectra. The regulation parameters and lock phase can be adjusted using an AU-program (Bruker C-program language) with the command xau set_loopval (after lock). The lists after running the AU-Prog give the option to select how strong the lock reacts on changes: lock.1: soft to lock.12 strong.

	1	2	3	4	5	6	7	8	9	10	11	12
Lock	119.3	115.4	110.2	107.2	104.1	99.7	96	92.6	89.6	86.0	83.9	82.2
gain												
Loop	-17.9	-14.3	-9.4	-6.6	-3.7	0.3	3.9	7.1	9.9	13.2	15.2	16.8
gain												
Loop	0.681	0.589	0.464	0.384	0.306	0.220	0.158	0.111	0.083	0.059	0.047	0.041
time												
Loop	20	30	50	70	100	160	250	400	600	1000	1500	2000
filter												

Table 1: Macros lock.x to set LOCK regulation parameters. Data from BSMS manual.

The command lock sets all the LOCK parameters and locks automatically on the selected solvent. The lock signal can be visualized in lockdisp. All LOCK parameters are listed in edlock.



The room temperature shims have to be optimized on every new sample to get optimal homogeneity and small line width of the NMR signal. The Bruker tool Topshim helps to adjust the shims, starting command: topshim gui. More information and tips for using Topshim:

- 1D Topshim optimizes only all the z-shims. The 1D Topshim can be used for all solvent. The program switches automatically to ¹H or ²H gradient shimming depending on the selected solvent during the lock procedure.
- 3D Topshim adjusts all the shims and is only for water samples.
- for shimming of x/y shims in topshim 1D use the tune option
- Z6 works on ¹H (less on ²H), check the water line in zgpr
- select parameters
 - o for Shigemi tubes: shigemi
 - o shimming Z7 & Z8: ordmax = 8
- 3D Topshim adjusts all the shims and is only for water samples.
 - o parameter: ordmax = 5
 - use before and after 1D topshim, after manual adjustment of x + y shim

Determine acquisition parameters **Pulse calibration** AUTOMATIC

The automatic ¹H pulse calibration using a single scan stroboscopic nutation experiment (P.S.C. Wu & G. Otting, Magn. Reson. 176, 115-119 (2005), description: <u>http://u-of-o-nmr-facility.blogspot.de/2010/03/fast-90-degree-pulse-determination.html</u>) is implemented in Topspin as an AU-Program (Bruker c-program language): pulsecal (for a water sample) or pulsecal sn (for all the other solvents).

MANUAL

The quality of the automatically ¹H pulse calibration depends on the spectrometer hardware and software. For that reason, the pulse calibration should also check manual following the procedure below:

•	read parameter set	\rightarrow	rpar PROTON
•	change pulse program to zg or zgpr	\rightarrow	pulprog zg
•	set the acquisition time to 100ms	\rightarrow	aq 100m
•	set number of scans to 1	\rightarrow	ns 1
•	set number of dummy scans to 0	\rightarrow	ds 0
•	set relaxation time delay to 1s	\rightarrow	d1 1
•	set receiver gain to 1	\rightarrow	rg 1
•	set proton pulse first short to 0.1 μ s	\rightarrow	p1 0.1
•	record a 1D spectrum	\rightarrow	zg
•	process the spectrum	\rightarrow	efp



- phase correction with only 0 order → apk0
- select the region around the solvent peak or other strong signal → dpl
- Optimization of the 360° pulse using the macro popt (parameter optimization), fill up the new window:

OPTIMIZE	GROUP	PARM.	OPTIMUM	START.	END.	NEXP	VARMOD	INC
Step by step	0	P1	ZERO	а	b	11	LIN	0.2

a = 4 times the pulse length from macro pulsecal, **b** = automatic

The result is stored in procno 999. The spectrum should show a zero crossing and the 360° pulse length is written in the title. If the zero crossing in the spectra is not obtained, change the popt parameter. Remember the ¹H-90° pulse length for the next experiments.

Prosol table

The prosol table is a list of standard pulse length and power level related to the Bruker pulseprogram library. The pulses are stored for each probe head, nucleus and channel. The prosol table can be opened with the command edprosol: new entries or change of entries can be done.

	w <u>H</u> eib							
		Saved Observ	e and Saved E	Decou	ple Prosol Parame	ter Set for:		
robe: Prodig	gy Z129773	3_0008 CPP TCI 60	00S3 H&F-C/I	N-D-(05 Z Select		Solvent: ge	neric 💌
		Observ	e		Deco	uple		
		1H	-	Nucle	us 13C	-		
		Observ	e		Deco	uple		
bserve Comm	ent: Defau	lt 1H obs 600		D	ecouple Comment	: Default 130	dec 600	
0 deg. Pulses	Square Pu	ulses Shape Puls	es Others	1				
	de contra de	Observe				Decouple		
	Nucleus	Pulse Width[µs]	Att. Lvl.[dB]	Set	Pulse Width[µs]	Att. Lvl.[dB]	Set Nucleus	
	1H	8.50	-10.00	0	8.00	-10.00	S 1H	
	2H	97.00	-13.80	0	97.00	-13.80	S 2H	
	13C	12.40	-22.55	0	12.40	-22.55	S 13C	
	15N	35.00	-24.77	0	35.00	-24.77	S 15N	
	19F	10.00	-10.41	0	10.00	-10.41	🚫 19F	
	Nucleus	Pulse Width[µs]	Att. Lvl.[dB]	Set	Pulse Width[µs]	Att. Lvl.[dB]	Set Nucleus	

The translation between prosol table and the pulse sequence definition of pulses, power levels and delay is defined in the relation files. The definition is indicated on the top of the pulse sequence by the command: prosol relations=<name>, e.g. prosol relations=<triple>. If in the pulse sequence no file is specified, the program uses the default prosol relation file, compare



file:///...topspin3.5home-directory/prog/docu/english/xwinproc/html_pp/792436363.html. The relation

files are stored in: ...TopSpin3.5home\conf\instr\spect\prosol\pulseassign. The command getprosol changes all pulses in the current data set to the values stored in the prosol table. In addition, the prosol table allows to recalculate the complete set of pulses for one nucleus,

e.g. if the ¹H is different from the standard pulse length, put a new value in and press the calc. button. Please do not store the changes. Only the spectrometer admin should do this. A better option is to use the command line options: getprosol Nuc P90 PL90 or write the new value in a macro, e.g. edmac ubi.

New: parameter can also copy to other parameter sets by command copypars.

Commands for the acquisition

- eda list all acquisition parameters
- ased only parameters listed which are necessary for using the actual pulse sequence
- rga automatic receiver gain adjustment
- zg delete the old data and starts the experiment
- go starts the experiment and add the data to the actual data set (works also in 2D)
- stop stops the experiment; losing the 1D data
- halt stops the experiment; store the 1D data on disk
- tr transfer data from acquisition memory to the disk
- gs pulsing of the pulse sequence without recording the data; for parameter optimisation

Commands for the 1D processing

- edp setup processing commands
- ft fourier transformation
- apk automatic phase correction
- apks automatic phase correction, faster
- abs
 base line correction and automatic integration
- sref
 referencing to solvent (solvent in eda)
- si number of points for processing

si = td * 0,5

- absf1, absf2 high field and low field values for abs in [ppm]
- absg, absl polynome used for abs, normal =5 and 3
- 2s nc_proc internal scaling factor

A selection of useful standard directories in TopSpin for

- pulse program: /<TopSpin_home>/exp/stan/nmr/lists/pp
- shape pulses: /<TopSpin_home>/exp/stan/nmr/lists/wave
- au-program: /<TopSpin_home>/exp/stan/nmr/au/src
- Parameter set: /<TopSpin home>/exp/stan/nmr/par
- User file were stored in additional user directory, e.g.
 - /<TopSpin_home>/exp/stan/nmr/lists/pp/user



Small molecules

Start a 1D 1H Experiment

¹H 1D-NMR Experiment

For setting up a standard ¹H 1D-NMR experiment:

rpar PROTON all & do all steps described in the general part. A useful test sample for the small molecules is 100mg Cholesteryl Acetat in C_6D_6 .

The PROTON parameter set use the pulse sequence zg30, where the 90° pulse is multiplied with 0.33. The smaller flip angle for the excitation pulse allows shorter relaxation delay (d1). Typical relaxation delays range from 1-5 times T₁. For best sensitivity choose the 'Ernst angle':

$$\cos \beta_{opt} = \exp(-TR/T_1)$$

(TR = pulse repetition time, TR = d1 + AQ; β_{opt} : optimum flip angle)

The consequences of too short relaxation delay are loss in sensitivity and artifacts in spectra. The pulse sequence with 30 degree pulses, e.g. zg30 minimized those problems. The pulse program code:

;zg30		
;1D sequence		
#include <avance.incl></avance.incl>	load definition file for NMR experiments	
1 ze	reset all memory on CPU	
2 30m	delay 30m	
d1	relaxation delay	
p1*0.33 ph1	proton 90° pulse multiplied with 0.33 at p	hase ph1
go=2 ph31	acquisition period	
30m mc #0 to 2 F0(zd)		0 is x axis
exit		1 is y axis
ph1=0 2 2 0 1 3 3 1	phase cycle of pulse p1	2 is -x axis
ph31=0 2 2 0 1 3 3 1	phase cycle of receiver	3 is -y axis
p1*0.33 ph1 go=2 ph31 30m mc #0 to 2 F0(zd) exit ph1=0 2 2 0 1 3 3 1 ph31=0 2 2 0 1 3 3 1	proton 90° pulse multiplied with 0.33 at p acquisition period phase cycle of pulse p1 phase cycle of receiver	hase ph1 0 is x axis 1 is y axis 2 is -x axis 3 is -y axis

All 1D/2D NMR pulse sequences are summarized in the TopSpin manual:

'Pulse Program Catalogue, 1D/2D'

Optimal Parameter for ¹ H-1D-I	NMR Experiment	
Acquisition time	aq	4s
Relaxation delay	d1	0.1
Spectral width	SW	20.0 ppm
Excitation frequency	o1 p	7.0 ppm
Receiver gain	rg	use rga
Number of scans	ns	32
Number of dummy scans	ds	8



Through choosing the combination of acquisition time and spectral width, the time domain points TD will be automatically calculated in TopSpin, because the parameters depend on each:

AQ = TD * DW= TD/2SWH (DW = dwell time, SWH = spectral width in Hz) For the 1D record a spectrum with large spectral width SW due the digital filtering will be neglected signals outside of spectral window. The combination of spectral width SW and excitation frequency o1p should be optimized for setting up a 2D-NMR spectra. The selected number of scans depends on the sample concentration.

¹³C 1D-NMR and DEPT Experiment

The natural abundance of the NMR active carbon isotope ¹³C is only 1 %. Therefore, is the sensitivity the main issue for ¹³C. In ¹³C-1D is the decoupling of ¹H necessary and results into two feature: simplify the multiplet structure and increase sensitivity due to Nuclear Overhauser Effect 'NOE'. The NOE depends also on the gyomagnetic ratio and can either be either negative or positive. For ¹³C is the NOE positive and decoupling also during the relaxation delay is useful. The different possibilities are summarized in table 3.

	Pulse program	Properties	Performance
$I \xrightarrow{d1} \qquad \qquad$	zg zg30	No decoupling	No signal enhancement, coupled spectrum
$I \xrightarrow{30^{\circ}} I \xrightarrow{1}$	zggd zggd30	,Gate Decoupling': decoupling during relaxation delay	Signal enhancement due to NOE effect, coupled spectrum
S CPD pl12			



$I \xrightarrow{d1} \downarrow \downarrow$	zgig zgig30	,Inverse Gate Decoupling': decoupling during acquisition	For quantitative analysis of ¹³ C- signals, decoupled spectrum
zgpg30 pli3-pli2 zgdc30 30° I d1	zgdc zgdc30 zgpg zppg30	Decoupling during entire experiment PG: ,Power Gate Decoupling' lower power during relaxation delay, pl13	For maximal signal intensity, decoupled spectrum

Table 3: pulse sequences for ¹³C 1D-NMR, pp- figures from Bruker manual, 'Pulse Program Catalogue, 1D/2D'.

Decoupling statements in the pulse sequence:

• •	
cw:	continuous wave decoupling
hd:	homo decouling
cpds1, cpds8:	composite pulse decoupling with CPD sequence 1, 8, synchronous mode
cpd1, cpd8:	composite pulse decoupling with CPD sequence 1, 8, asynchronous mode
do	switch decoupling off

Example: part of zgpg (F1-Channel: ¹³C & F2-Channel: ¹H)

•••••		
10u pl13:f2	change power level on chan	nel 2, soft decoupling power
d11 cpd2:f2	start decoupling	
DELTA		
4u do:f2	switch decoupling off	
10u pl12:f2	change power level on chan	nel 2, strong decoupling power
100m cpd2:f2	start decoupling program cp	d2
p1*0.33 ph1		
go=2 ph31		
30m do:f2 pl13:f2 mc #0 to	2 F0(zd)	switch decoupling off
exit		

The most popular pulse sequence for carbon 1D is zgpg30, due maximal signal intensity. Do rpar C13CPD all (pp=zgpg30) and getprosol. Further optimization like pulse calibration are not necessary for hetero nucleus. Use maximum receiver gain and start the experiment.



Due the decoupling in a ¹³C 1D-NMR the information how many protons bound to the carbons is neglected. Polarization transfer technique like DETP bring back this information. The DEPT starts on ¹H and transfer the magnetization to ¹³C, which increase the sensitivity compared to the standard ¹³C 1D-NMR. DEPT Experiment can easily start by:

rpar C13DEPT135 and getprosol (or better getprosol 1H P90 PL90, if you have calibrated the ¹H-pulse). The parameter set based on the pulse sequence deptsp135, which use a adiabatic composite pulse for refocusing. To compare how the sensitivity can profit from using adiabatic shape in pulse sequence generate a new data set and change the pulse sequence to dept135 and re- measure the DEPT.

Standard 2D Experiments

A wide range of possible and useful pulse sequences for running a 2D COSY experiment are included in TopSpin. The important NMR experiments for analyzing small molecules is summarized in the data set for CMSse which use optimized pulse sequences. The command rpar CMCse* show all experiments

🍬 Parameter	Sets: rpar CMCse*			×
File Options	<u>H</u> elp	Source = C	:\Bruker\TopSpin3.5\	exp\stan\nmr\par ~
Find file names	✓ CMCse*	Exclude:	Clear	
Class = Any	→ Dim = Any →	Show Recommended		
Type = Any	✓ SubType = Any ✓	SubTypeB = Any ~	Reset Filters	
CMCse_13C	CMCse_15NHMBCf2	CMCse_15NHSQCf2	CMCse_1H	CMCse_ADEQ
CMCse COSY	CMCse H2BC	CMCse HMBC	CMCse HSQC	CMCse INAD
				Read Close

CMCse is an extra Bruker program for automatically analyzed NMR spectra of small molecules which need an additional license)

Or select the recommended data sets:

🖕 Parameter Sets: rpar				E ×	
<u>F</u> ile <u>O</u> ptions <u>H</u> elp			Source = C:\Bruk	xer\TopSpin3.5pl2\exp\stan\nmr\par ~	
Find file names V enter any st	ring, *, ? Exclude:	Clear			
Class = Any · Dim = Any	y 🗸 Show Recommended				
Type = Any V SubType	e = Any \sim SubTypeB = Any \sim	 Reset Filters 			
C13CPD	C13DEPT135	COSYGPDFPHSW	COSYGPSW	HMBCETGPL3ND	
HMBCGP	HMBCGP_15N	HSQC_TOCSY	HSQC_TOCSY_ADIA	HSQCEDETGPSISP	
HSQCEDETGPSISP_ADIA	HSQCETGP_15N	HSQCETGPSISP	HSQCETGPSISP_ADI	A MLEVPHPR	
MLEVPHSW	NOESYPHPR	NOESYPHSW	PROTON	ROESYPHPR	
ROESYPHSW	WATERSUP				
				Read Close	



2D COSY

rpar CMSseCOSY

getprosol 1H P90 PL90

The parameter set based on pulse sequence cosygpmfppqf. A pulse sequence using gradient pulses for selection with multiple quantum filter according to gradient ratio.

The sample depending parameters are acquisition time (AQ), spectral widths (SW) and excitation frequency (o1p) should set in eda. For the indirect dimension the parameter SW1 = SW. The minimum value for TD1, number of points in the indirect dimension should in minimum set to 384, because the polarization transfer in COSY is happening during the t1 evolution period. That is also the reason why the receiver gain adjustment where doing in the ¹H 1D-NMR. Due the gradient selection with multiple quantum filter the experiment works only with 1 scan.

For suppression the zero quantum coherence by the z-filter (M.J. Thrippleton & J. Keeler, Angew. Chem. Int. Ed. 42, 3938-3941 (2003)) change:

- pulprog cosygpphzfzs
- 1 FNMODE States-TPPI
- gpz0 10
- NS 4

selective 1D-NOESY experiment

Introduction

This experiment consist of three parts:

- Selective excitation of the selected resonance using the SPFGE block.
- Mixing period consisting of the basic 90⁰(1H)-delay-90⁰(1H) block in phase polarization transfer to other spins via NOE. Purging gradients are usually applied during the mixing period in order to remove any residual transverse magnetization.
- Proton detection as usual.

Figure 5.21



Steps for setup a selective NOESY:

• record a ¹H 1D-Experiment



- select separated signals and determine the chemical shifts
- generate a new dataset
- rpar SELNOGP all
- getprosol 1H P90 PL90
- **o1p** = set to the chemical shift of the selected signal of step 2
- d8 = 300 ms or better the T₁-time of the selective signal of step 2

For each selected signal a new EXPERIMENT must have recorded. Be aware that the sign between selected signals and the NOE-signals is the opposite. Instead of using the standard pulse sequence selnogp, use selnogpzs with better artefact suppression.

selective 1D-TOCSY experiment

- This experiment consist of three parts:
- · Selective excitation of the selected resonance using the SPFGE block.
- Mixing period to achieve in phase polarization transfer to other spins. This is usually achieved by applying some isotropic mixing sequence like MLEV, WALTZ or DIPSI pulse trains. This in-phase transfer avoids possible cancellation when the coupling is poorly resolved.
- Proton detection as usual.

Figure 5.38



Steps for setup a selective TOCSY:

- record a ¹H 1D-Experiment
- select separated signals and determine the chemical shifts
- generate a new dataset
- rpar SELMLGP all and change pp to seldigp
- getprosol 1H P90 PL90
- **o1p** = set to the chemical shift of the selected signal of step 2
- d9 = 20ms for only one transfer step
- or set d9 = 80ms for more transfer step

For each selected signal a new EXPERIMENT must have recorded. The selected signal and the TOCSY-signals have the same sign. Instead of using the standard pulse sequence seldigp, use seldigpzs with better artefact suppression.

2D HSQC / H2BC



Compare the different HC-correlations:

Experiment A Experiment B Experiment C rpar HSQCETGPSI rpar CMCse_HSQC rpar C1MCseH2BC getprosol 1H P90 PL90 AQ, o1p, SW: optimized value from the 1D o2p 75ppm 1 SW 150ppm 1TD 64 NS 4

Compared the results, which HC correlation is the most useful experiment.

2D HMBC / sel HMBC Set up: rpar CMCse HMBC getprosol 1H P90 PL90

The pulse sequence hmbcetgpl3nd uses three-fold low-pass J-filter to suppress one-bond correlations (D.O. Cicero, G. Barbato & R. Bazzo, J. Magn. Reson. 148, 209-213 (2001).

AQ, o1p, SW: optimized value from the 1D o2p 100ppm 1 SW 250ppm 1TD 64 NS 8

The processing command for the HMBC is xfb and xf2m. For running a selective HMBC change the pulse sequence to shmbcctetgpl2nd and change the recommended gradient values. Start the shape tool (stdisp) in the ¹³C 1D-data set and select a region of interest. Use the button in shape tool to find the right pulse length for the shape

🖕 Edit regions	: for Dataset: <2016-02-2	6-G-NMR-C 🖪 🗙
Use advanced	Imode	
🗹 Use same sha	pe for all regions	
Carrier frequency	O1 = 15091.29 Hz	
Left Rig	Shape	Type of rotation
68.25 12.82	Gauss Cascade Q3 $$	Inversion (Iz->-Iz) \sim
	<u>O</u> K	<u>Apply</u> <u>Cancel</u>

Set the pulse length (p32) and offset (spoffs32) of the shape. In a second experiment set the offset of shape (spoffs32) to 0 and change the middle of the 13C-frequency with smaller spectral widths in F1. Compare all three HMBC spectra.



Bio-NMR

Start a 1D 1H Experiment

WATERGATE

rpar ZGGPWG all

getprosol 1H P90 PL90

Pulse sequence code: z	ggpwg	
prosol relations= <triple< td=""><td>></td><td>7000100</td></triple<>	>	7000100
<pre>#include <avance.incl></avance.incl></pre>		zggpwg
#include <grad.incl></grad.incl>		-x
"p2=p1*2" "d12=20u"		
1 ze		sp1
2 30m		G _z
d1		G1 G1
10u pl1:f1		
p1 ph1		
50u UNBLKGRAD		
p16:gp1	first gradient dephases all the coher	ences, gpnam1=SINE.100, gpz1=20%,
p16=1ms		
d16 pl0:f1	recovery delay, set power for F1-Chann	el to pl0
(p11:sp1 ph2:r):f1	sel 90° pulse on water, p11 = 1ms, spnar	m1=Sinc1.1000, sp1 calculated from p1
4u		
d12 pl1:f1 (p2 ph3)	set power back to hard power of proto	n pulse
4u		
d12 pl0:f1		
(p11:sp1 ph2:r):f1	second selective 90° pulse for v	water, the phase 2 can be optimized in
gs-mode		
46u		
p16:	gp1 second gradient re-phases all the i	esonances exceptto the on resonance
solvent resonance		
d16 4u BLKGRAD go=2 ph31 30m mc #0 to 2 F0(zd)		
EXIT	(M. Pio	tto, v. Saudek, v. Sklenar, J.Biomol.,
1992, 2, 661)		
ph1=0 2 ph2=0 0 1 1 2 2 3 3 ph3=2 2 3 3 0 0 1 1 ph31=0 2 2 0		
p101-0 2 2 0		



Excitation Sculpting

rpar ZGESGP all

getprosol 1H P90 PL90

Pulse sequence code: zgesgp prosol relations=<triple>

#include <Avance.incl>

#include <Delay.incl>

"p2=p1*2" "d12=20u" "TAU=de+p1*2/3.1416+50u" "acqt0=0" baseopt_echo 1 ze 2 30m d12 pl1:f1 BLKGRAD d1 p1 ph1 50u UNBLKGRAD p16:gp1 d16 pl0:f1 (p12:sp1 ph2:r):f1 calculated from p1 4u d12 pl1:f1 p2 ph3 4u p16:gp1 refocussing gradient d16 TAU p16:gp2 d16 pl0:f1 (p12:sp1 ph4:r):f1 4u d12 pl1:f1 p2 ph5 non selective 180° pulse 4u p16:gp2 refocussing gradient d16 go=2 ph31 30m mc #0 to 2 F0(zd) 4u BLKGRAD exit ph1=0 ph2=0 1 ph3=2 3 ph4=0011 ph5=2 2 3 3 ph31=0 2 2 0





Start first gradient echo, p16=1ms, gpnam1=SINE.100, gpz1=31% recovery delay, set power for F1 to pl10=120Db selective 180° pulse for water, spnam1=Sinc1.100, p12=2ms, sp1 set power F1 channel back to pl1

non selective 180° pulse

start second gradient echo, p16=1ms, gpnam2=SINE.100, gpz2=11%

selective 180° pulse for water

(T.-L. Wang, A. J. Shaka, J.Mag.Res., 1995, A112)



More Water suppression sequences:

- zggpw5,
- zggpjrse, Jump-Return Echo, the best pulse sequence for analysis of RNA.
- zgesgppe, with perfect Echo, reduces the artifacts in 1D.

set up 2D ¹⁵N-¹H correlation, 600 MHz

Setting up a 2D ¹⁵N-¹H correlation NMR experiments, like HSQC, starts with loading the Bruker

standard parameter set using the commands:

rpar <<PARAMETER SET>> getprosol 1H P90 PL90

For ¹⁵N- and ¹³C-labeled Protein the ZGOPTNS –DLABEL_CN should be activated to enable the ¹³C decouplingduring the ω_1 -evolution period.

General starting parameters for a protein-project at 600 MHz are:

Aq (¹ H) = 71r	ns	Aq (¹⁵ N) = 105ms
SW (¹ H) = 12	opm	SW (¹⁵ N) = 40ppm
o1p (¹ H) = 4.	7ppm	o3p (¹⁵ N) = 117ppm
d (d26) = 2.6	31ms (cnst4 = 95 Hz)	
¹⁵ N dec.: pcp	d3 = 220µs	CPD_prog = <i>grarp4</i>
o2p (¹³ C) = 101ppm	¹³ C dec.: p14 = 500µs	adiabatic pulse = Crp60,0.5.20.1

These general parameters should be optimized for each protein-project and remembered for the 3D experiments below.

in dualition, ener	en une parameters se			
	HSQC	HMQC	TROSY	TROSY (sel. ¹⁵ N)
Parameter set	HSQCFPF3GPPHWG	SFHMQCF3GPPH	B_TROSYETF3GPSI	B_TROSYETF3GPSI
Pulse program	hsqcfpf3gpphwg	sfhmqcf3gpph	b_trosyetf3gpsi.2	b_trosyetf3gpsi.3
DS	32	32	32	32
NS	8	32	32	32
d1 [s]	1	0.2	0.2	0.2
¹ H selective pulses	sp1 (90°pulse) Sinc1.1000 p11 → 1 ms	sp23 (120°pulse) Pc9_4_120.1000 p39 \rightarrow 2.4 ms	sp26 (180°pulse) Repurp.1000 p42 → 1.4 ms	sp25 (90°pulse) Pc9_4_90.1000 p41 → 2.2 ms
		sp24 (180°pulse) Rsnob.1000 p40 → 0.8 ms	sp28 (90°pulse) Eburp2.1000 p43 → 1.7 ms sp29 (90°pulse) Eburp2tr.1000 p43 → 1.7 ms	sp26 (180°pulse) Repurp.1000 p42 \rightarrow 1.4 ms sp28 (90°pulse) Eburp2.1000 p43 \rightarrow 1.7 ms sp29 (90°pulse) Eburp2tr 1000
				p43 → 1.7 ms
¹⁵ N selective pulses				sp39 (180°pulse) Reburp.1000 p56 → 1.6 ms

In addition, check the parameters below:



sp40 (180° pulse) Bip720,50,20.1 p56 \rightarrow 0.5 ms^A

Table2: Parameters for 2D $^{15}\text{N-}^{1}\text{H}$ correlation experiments. ^ 500 μs at all field strength

New option in Topspin 3.5: Setting ZGOPTNS to -DCALC_SP (or -DCALC_SP -DLABEL_CN for ¹³C

& ¹⁵N labeled proteins) allows to automatically calculate all the band-selective proton pulses, based

on cnst54 and cnst55 and the manually determined p1

cnst54: H(N) chemical shift (offset, in ppm), Protein: 8.5 RNA: 12.3 cnst55: H(N) bandwidth (in ppm), for both 4.8

Shape pulses

Calculation of radiofrequency field strength γB_1

 The magnetization vector precesses around the radiofrequency field B₁ according to

$$\omega = 2 \pi v = \gamma B_1 \tag{1}$$

- The angle precession q (in rad), is proportional to the pulse width p: $\theta = \gamma B_1 * \pi \tag{2}$
- For a 90° pulse ($\theta = \pi/2$) and Equation (1):

$$\gamma B_1 = 1/(4 * \pi)$$
 (3)

Types of shape pulses

- 90° pulses (excitation, $I_z \rightarrow I_y$)
- 180° pulses
 - Inversion pulses, $I_z \rightarrow -I_z$
 - Refocusing pulses, $I_y \rightarrow -I_y$
- Very narrow excitation
 - Solvent suppression (e.g. Sinc, Gauss, Square pulses)
- Excitation / Inversion of a region Ca, C'...
 - o Gaussian Cascades
 - 180° pulses: Q3, G3 90° pulses: Q5, Q5tr, G4, G4tr
- Broadband Inversion/ Refocusing
 - Inversion, $I_z \rightarrow -I_z$
 - Adiabatic pulses (e.g. smoothed Chirp pulses)
 - Broadband Inversion pulses (BIPs)
 - Refocusing, $I_x \rightarrow -I_x$ or $I_y \rightarrow -I_y$
 - Composite adiabatic pulses (e.g. Composite Chirp)

Nomenclature for shapes in Topspin:

- Standard shapes, e.g. Q5tr.1000 →
 - o Shape (Q5 time reversed)
 - Number of points (1000)
- Adiabatic shapes, e.g. Crp60, 05,20.1



- Shape and sweep width (Chirp, 60 kHz)
- o duration (0.5 ms)
- Smoothing and number of points (20% and 1000)

Instead of using band-selective pulses, for protein NMR selective rectangular pulses can be applied due the large chemical shift difference between CO and CA. The excitation profile with maximum on CA and minimum on CO is calculated as follows:

90° pulse:

$$\beta^{\text{eff}} = 4*\beta = 360^{\circ} \Longrightarrow (\omega_1^{\text{eff}})^2 = (4\omega_1)^2 = \omega_1^2 + (\Delta\Omega)^2 \Longrightarrow \omega_1 = \Delta\Omega/\sqrt{15}$$

$$\tau_{\mathrm{p}(90^\circ)}=\beta/\omega_1=\pi/2*\sqrt{15}\,/\,(2\pi*\Delta\Omega)$$

180° pulse:

$$\beta^{\text{eff}} = 2*\beta = 360^{\circ} \implies (\omega_1^{\text{eff}})^2 = (2\omega_1)^2 = \omega_1^2 + (\Delta\Omega)^2 \implies \omega_1 = \Delta\Omega/\sqrt{3}$$

$$\tau_{\mathsf{p}(180^\circ)} = \beta/\omega_1 = \pi * \sqrt{3} / (2\pi * \Delta \Omega)$$

Both equations can be implemented in a pulse sequence using cnst21 (CO chemical shift) and cnst22 (CA chemical shift) with the following definitions:

90°-pulse: "pxx=3.873*1000000/((cnst21-cnst22)*bf2*4)" 180°-pulse "pxx=1.732*1000000/((cnst21-cnst22)*bf2*2)"

Setup 3D Experiment

First read a standard parameter set, e.g. for the protein backbone 3D experiments

rpar <<PARAMETER SET>> getprosol 1H P90 PL90

set the optimized general parameters from the HSQC

- protein backbone 3D experiments based on WATERGATE: HNCOGPWG3D, HNCACOGPWG3D, HNCAGPWG3D, HNCOCAGPWG3D, HNCACBGPWG3D, HNCOCACBGPWG3D
- protein backbone 3D experiments based on Best-method:
 B_HNCOGP3D, B_HNCACOGP3D, B_HNCAGP3D, B_HNCOCAGP3D;
 B_HNCACBGP3D, B_HNCOCACBGP3D
- protein backbone 3D experiments based on Best-TROSY-Method:
 B_TRHNCOGP3D, B_TRHNCACOGP3D, B_TRHNCAGP3D, B_TRHNCOCAGP3D;
 B_TRHNCACBGP3D, B_TRHNCOCACBGP3D

Which 3D methods is the best option for protein assignment depends on the size of the molecule. This can be experimentally tested by recording the carbon plane in a 3D HNCO. This



can be done by setting the number of points in the evolution period of ^{15}N (F2) to 1 and ^{13}C (F3) to 128, td.

🍬 TD		B×
Size of fid (F3, F2, F1)		
TD 1024	11	128
		<u>O</u> K <u>C</u> ancel

Processing the 2D Plane in a 3D dataset.

• Checking reverse statement during transformation: reverse

🍬 REVER	RSE		E ×
Reverse sp	pectrum during	transform (F3, F2, F1)	
REVERSE	FALSE	∽ FALSE	TRUE
			<u>O</u> K <u>C</u> ancel

Depends on the pulse sequence, hncogpwg3d: FALSE; b_hncogp3d.2 & b_trhncogp3d.2:

TRUE

- 2D processing command FT: xfb
 - Enter a new PROCNO for 2D data: e.g. 13 (store the carbon-plane in procno 13 and ¹⁵N plane in procno 23).
- 2D processing command for phase correction: apk2d
- Compare the 2D spectra:
 - read out the 1D projection, by rhpp 2 (read horizontal positive projection and store it in procno 2)

Additional Information

All implemented NMR-experiments in TopSpin are described in the manual 'pulse program

catalogue BIO'.

In TopSpin a wide range of NMR experiments are implemented as parameter set which can

easily use. A selection of important experiments follows below:

Parameter sets for protein experiment with best method: rpar B_*

🍦 Parameter Sets: rpar B_*					×
<u>File</u> Options <u>H</u> elp			Source = C:\Bruker\To	pSpin3.5\exp\stan\nmr\par	\sim
Find file names $\ \ > \ B_*$	Exclude:	Clear			
Class = Any V Dim = /	Any 🗸 🗌 Show Recommen	ided			
Type = Any V SubTy	rpe = Any \checkmark SubTypeB = A	Any V Reset Filters			
B HNCACBGP3D	B_HNCACBIGP3D	B_HNCACOGP3D	B_HNCACOGP4D	B_HNCAGP3D	
B_HNCAIGP3D	B_HNCOCACBGP3D	B_HNCOCACBGP4D	B_HNCOCAGP3D	B_HNCOCAGP4D	
B_HNCOGP3D	B_HNCOIGP3D	B_HSQCETF3GPSI	B_TRHNCACBGP3D	B_TRHNCACBIGP3D	
B_TRHNCACOGP3D	B_TRHNCAGP3D	B_TRHNCAIGP3D	B_TRHNCOCACBGP3D	B_TRHNCOCAGP3D	
B_TRHNCOGP3D	B_TRHNCOIGP3D	B_TROSYETF3GPSI	B_TROSYF3GPPH		
				Read CI	ose



Parameter sets for 4 dimensional NMR experiments: rpar *4D*

🍬 Parameter Sets: rpar *4D	×				×
<u>F</u> ile <u>O</u> ptions <u>H</u> elp			Source = C	C:\Bruker\TopSp	in3.5\exp\stan\nmr\par ~
Find file names V 4D*	Exclude:	Clear			
Class = Any \checkmark Dim = A	ny 🗸 🗌 Show Recommende	ed			
Type = Any V SubTy	pe = Any \sim SubTypeB = Any	y V Reset Filters			
B HNCACOGP4D	B_HNCOCACBGP4D	B_HNCOCAGP4D	CBCACONHGP	WG4D	HBHACBCANHGPWG4D
HBHACONHGPWG4D	HNCACOGP2H4D	HNCACOGP4D	HNCACOGPWG	2H4D	HNCACOGPWG4D
HNCOCAGP2H4D	HNCOCAGP4D	HNCOCAGPWG2H4D	HNCOCAGPWG	64D	HSQCNOESYHSQCCCGP4D
HSQCNOESYHSQCCNGP4D	HSQCNOESYHSQCNCGP4D	HSQCNOESYHSQCNNGP4D	TRHNCOCAET	GP2H4D	TRHNCOCAETGP4D
TRHNCOCAGP2H4D	TRHNCOCAGP4D				
					Read <u>C</u> lose

Parameter sets for ¹³C direct detection: rpar C_*

🖕 Parameter Sets: rpar C_*					×
<u>F</u> ile <u>O</u> ptions <u>H</u> elp			Source = C:\E	Bruker\TopSpin3.5\exp\stan\nmr\par	\sim
Find file names $\ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \ \$	Exclude:	Clear			
Class = Any V Dim = Ar	ny V Show Recommended	t			
Type = Any V SubType	e = Any \checkmark SubTypeB = Any	✓ Reset Filters			
C CACO	C_CACO_IA	C_CACO_S3	C_CAN_IASQ	C_CAN_MQ	
C CAN MQ.2	C CANCO IA3D	C CANCO IA3D.2	C CANCOI IA3D	C CBCACO IA3D	
C_CBCACO_S33D	C_CBCACON_IA3D	C_CBCANCO_IA3D	C_CCCO_IA3D	C_CCCO_S33D	
C_CCCON_IA3D	C_CCFLOPSY16	C_CCFLOPSY16_CT	C_CCFLOPSY16_CT	IA C_CCFLOPSY16_IA	
C_CCNOESY	C_CCNOESY_CT	C_CCNOESY2	C_COCA	C_COCA_IA	
C_COCA_MQ	C_COCA_MQ.2	C_CON_IASQ	C_CON_MQ	C_CON_MQIA	
C_CON_SQ	C_COSY	C_COSY_CT	C_COSY2_CT	C_HACACO_3D	
C_HCACO_IA3D	C_HCACO_S33D	C_HCAN_IA3D	C_HCANCO_IA3D	C_HCANCOI_IA3D	
C_HCBCA_IA3D	C_HCBCACO_IA3D	C_HCBCACO_S33D	C_HCBCAN_IA3D	C_HCCFLOPSY16_3D	
C_HNCA_IA3D	C_HNCACO_IA3D	C_HNCACO_S33D	C_HNCO_IA3D	C_HNCOCA_IA3D	
C_HNCOCA2_IA3D					
				Read <u>C</u>	lose

Parameter sets for nucleic acids: rpar NA_*

🍦 Parameter Sets: rpar NA_*								
<u>F</u> ile <u>O</u> ptions <u>H</u> elp			Source = C:\Bruker\TopSp	oin3.5\exp\stan\nmr\par ∨				
Find file names \sim NA_*	Exclude:	Clear						
Class = Any V Dim = Any V Show Recommended								
Type = Any V SubTy	pe = Any \checkmark SubTypeB = Ar	ny ~ Reset Filters						
NA C6NOESYHSQCGP3D	NA GHCCHFWDIGP3D	NA H56C56C4N3H 3D	NA H56C56C4N3H 3D2	NA H5C5C4N3H 3D				
NA H68C68N19C42	NA HCCHECGP3D	NA HCCHFWDIECGP3D	NA HCCHFWDIECGPJR3D	NA HCCHFWDIGP3D				
NA_HCCNHDIGPWG3D	NA_HCNCCHCOMQ3D	NA_HCNCCHDIMQ3D	NA_HCNCHGPJRPHSP	NA_HCNCHGPMQSP3D				
NA_HCNETGPSI3D	NA_HCNETGPSISP3D	NA_HCNHGPPH19	NA_HCNMQ3D	NA_HCNMQGPPHPR				
NA_HCNMQSP3D	NA_HCPDIETGP3D	NA_HCPDIETGPSI3D	NA_HCPDIGP3D	NA_HCPDIGPJR3D				
NA_HCPETGPSI3D	NA_HCPQETGPSI.1	NA_HCPQETGPSI.2	NA_HNC6C5ETGPSI	NA_HNCCH3D				
NA_HNNCOSYGPPHSPWG	NA_HNNCOSYGPPHWG	NA_HPCCOCTETF3GP3D	NA_HPCCOETF3GP3D	NA_HPCCTCO3D				
NA_HPCOSYPHPR	NA_HPCTCO	NA_HPCTCOJR	NA_HPDI	NA_HPDICO3D				
NA_HPDINO	NA_HPDINO19	NA_HPDINO193D	NA_HPDINO3D	NA_HSQCCTETGPSP.2				
NA_HSQCETF3GPXY	NA_HSQCETGPJCSP	NA_HSQCF3GPJRPHXY	NA_HSQCF3GPNOPHXY	NA_HSQCFPF3GPPHWG				
NA_NOESYGPPHPRXF	NA_NOESYGPPHXF19	NA_NOESYHSQCCTETGP3D	NA_NOESYHSQCF3GPWG3D	NA_PCCHCO				
NA_PCCHDI	NA_PCCHDI2	NA_PFIDSETGPSI	NA_TRHCCHCO3D	NA_TRHCNCHMQSP3D				
NA TRHCNCHMQSPWG3D	NA_TRHCNCHMQSPWG3D2	NA_TRHCNET	NA_TRHCNETSI3D	NA_TRHCNPH				
NA_TRHCNPH3D	NA_TRHNNCOSYGPPHSPWG	NA_XHCOCTETF3GP	NA_XHCOETF3GP	NA_ZGGPJRSE				
				Read Close				



NUS – Non Uniform Sampling

The general concept of NUS is summarized in TopSpin manual 'NUS Parameters'. The NUS reduced for a multidimensional NMR-experiment the number of points for the indirect dimension by using algorithm. For using NUS in setup NMR-Experiments only minor changes of the parameter setup are necessary. First, you set up the 2D, 3D or 4D NMR experiments in the described way before. The NUS will be activated for the data set by selecting the parameter FnTYPE as non-uniform_ sampling in eda.

\$ <mark>1</mark>	2016-02-23-G-N	MR-Ubiquitin7 12 1	I d:\data\ric\nmr					
5	Spectrum ProcPars AcquPars Title PulseProg Peaks Integrals Sample Structure Plot Fid							
M								
E	xperiment		F4	F3	F2	F1	Frequency axis	^
R	/idth eceiver	S Experiment						
N	ucleus	PULPROG	b_hncocacbgp4d	De_hncocacbgp4d E DQD v non-uniform sampling v			Current pulse program	
	urations	AQ_mod	DQD ~				Acquisition mode	
	rogram	FnTYPE	non-uniform_sampling				nD acquisition mode for 3D etc.	
P	robe	FnMODE	traditional(planes) full(points) non-nunform_sampling			Acquisition mode for 2D, 3D etc.		
L	ists	ProjAngle [degree]				Angle for projection-spectroscopy		
N	US	TD				Size of fid		
V	/obble	DS	T6	rojection-spectroscopy ხ				

In the category NUS of eda the number NusAMOUNT [%] selected how many total points are recorded for the indirect dimension.

1 2016-02-23-G-N	MR-Ubiquitin7 12 1	l d:\data\ric\nmr					
Spectrum ProcPa	cPars AcquPars Title PulseProg Peaks Integrals Sample Structure Plot Fid						
🔊 Л S 🕇 🖾	s ⊎ tŽ., ♥ C 🚜						
Experiment Width	NUS (Non Uniform Sampling) parameters					^	
Receiver	NUS Help					Show NUS help	
Nucleus	NusAMOUNT [%]	25				Amount of sparse sampling	
Durations	NusPOINTS	1536				Number of hypercomplex points in indirect dimension	
Program	NusJSP [Hz]		0	0	0	J-coupling	
Probe	NusT2 [sec]		1	1	1	T2 relaxation	
Lists	NusSEED	54321				Random generator seed	
NUS	NUSLIST	automatic				Name of loopcounter list for NUS (Non Uniform Samp	
Wobble		Calculate				Calculate point spread function	
Automation		Show				Display NUS point spread	

The choice of NusAMOUNT [%] compromised between as less as possible points and getting a multidimensional NMR spectrum without artefacts. From practical experiences is used as recommended values 50% for 2D NMR and 25% for 3D & 4D NMR. The parameter NusT2 [sec] weight the calculated points combination for the indirect dimension. If all NusT2 1 yields in equal distribution of all points. For **protein NMR experiment** the suggestion of NusT2:

¹H(aliphatic):0.02 ¹⁵N: 0.05 ¹³C: 0.5,

which weight the points distribution more to beginning of the experiment, shown by press Show button. Result:





Processing of NUS data.

TopSpin recognized automatic if the date were recorded with NUS and do the processing like usually. For processing of NUS data two algorithm are implemented in TopSpin: **MDD** and **CS**. Both need an additional license for processing, but not for recording the data. Which algorithm are used can define by the parameter Mdd_mod in edp, showing below:

Spectrum ProcPars AcquPars Title PulseProg Peaks Integrals Sample Structure Plot Fid								
n S 12. M ♥ ∰								
	Reference Window	NUS (Non Uniform Sampling) parameters						^
	Phase	Mdd_mod	mdd 🗸 🗸			MDD mode		
	Baseline	MddCEXP	mdd	FALSE ~	FALSE ~	RMDD/MDD flag		
	Fourier	MddCT_SP	CS	FALSE ~	FALSE ~	Constant time		
	NUS Reak	MddF180		FALSE ~	FALSE ~	Delayed sampling flag		
	Automation	MddNCOMP	0			Number of components		
	Miscellaneous	MddPHASE		0	0	Phase		
	User	MddSRSIZE [ppm]	0			Sub region size		

General, in the most cases the **CS** algorithm give a better result, but the processing takes longer. For that reason, optimized the other processing parameter with **MDD** and select **CS** for final processing. In TopSpin 3.5 the **CS** algorithm is included for free by processing 2D data.