

SPARKY

Graphical NMR assignment and integration program
for proteins, nucleic acids, and other polymers



Henry Jonker

Center for Biomolecular Magnetic Resonance (BMRZ)
Institute for Organic Chemistry and Chemical Biology
Goethe University - Frankfurt am Main (Germany)

NMR data analysis with Sparky (NMR Intensiv kurs 2016)

Sparky

❖ Developed & maintained by Donald Kneller and Tom Goddard
(University of California, San Francisco)

- freely available for all common operating systems (Windows, Linux, Unix and Macintosh)
- www.cgl.ucsf.edu/home/sparky/ (download & online manual)
- tech.groups.yahoo.com/group/nmr_sparky/ (discussion group)

❖ Works simultaneously with any number of 2D, 3D or 4D spectra

❖ Can use all commonly processed NMR spectra from

- XWinNMR, Topspin (Bruker), VNMR (Varian), NMRPipe, Felix
- easy to convert and to edit the properties of the spectra


❖ NMRFAM-Sparky

- development version with additional tools by and for NMRFAM
(University of Wisconsin, Madison)

Directory setup:

❖ Essential folders and files for sparky datasets

- data
 - processed NMR spectra in sparky (.ucsf) format
- save
 - spectra (.save) files with all the display settings, peaks, assignments, labels, ...
- lists
 - output for peak and assignment lists, spectra (postscript), ...
- projects
 - session project (.proj) files to load multiple spectra with their relation settings
- sparky.bat (windows)
 - not essential, but usefull executable to the start sparky project (copy from sparky/bin folder and add the variable: SPARKYHOME=.)
 - check the python version for different versions of sparky



Data conversion:

❖ Convert spectra to .ucsf data:

- example: **bruk2ucsf** 2rr n15hsqc.ucsf
 - will create a n15hsqc.ucsf file from the xwinnmr/topspin 2rr file
 - also: **vnmr2ucsf** for Varian spectra and **pipe2ucsf** for NMRPipe spectra

❖ Show spectrum information:

- example: **ucsfdata** hnca.ucsf
 - type only 'ucsfdata' to get a full list of all the options

❖ Edit spectrum settings (eg. axis names, frequency, offset,)

- example: **ucsfdata -a2** HN n15hsqc.ucsf
 - will rename the w2-axis to 'HN')
- example: **ucsfdata -p2 -o** n15noesy-2dhh.ucsf n15noesy.ucsf
 - makes a 2D ^1H -NOESY projection from the 3D ^1H ^1H ^{15}N -NOESY-HSQC dataset


Menu & Pointer

❖ Commands are available via the menu and 2-letter shortcuts

- right-click in spectrum for direct menu access

❖ The mouse pointer


- select function or directly change with keys F1-F12



Backup & Print

❖ Automated backups

- set the backup time period in the preferences dialog (pf)
- the .save file is copied to .save.BAK
- the .proj file is copied to .proj.BAK
- only one backup is kept




❖ Undo (eu) function

- works only one step back

❖ Print a spectrum (lt)

- select the spectrum and save as a postscript file
- print this file (use ao. ghostview)




View mode

❖ View options (vt) and contour levels (ct)

- change the look of the spectrum.
- peaks and the label sizes can be adapted (oz)


crosshair
(also connected to
other spectra views)



View mode

❖ Change aspect ratio

- widen or narrow the spectrum view



❖ Adapt the assignment format to change the peak labels

- %a1-%a2 full atom names and group residues G16N-HN
- %A1-%A2 force the groups to be displayed (noesy) G16N-G16HN
- %G1 show only residue name (of first axis) G16

Move, Zoom, Duplicate, Hide & Center

❖ Mouse pointer mode or shortcuts to zoom and center

- zoom in (zi), zoom out (zo)
- show full spectrum (zf), zoom to previous state (zp)
- center view (vc) on peak position

❖ Additional views of the same spectrum

- create a duplicate view (vd)
- useful to assign the other side of the diagonal in a 2D NOESY (rotate with xr)
- take care not to remove the original view instead of the duplicate!
(a bad situation will arise in which Sparky commands may not work)




❖ Manage spectrum views

- hide a spectrum view (vh)
- show the spectrum by selecting in the view list (pv) or in the window menu

Synchronization

❖ Synchronize (yt) the axes of 2 or more selected spectra

- the views of the connected 2D and/or 3D spectra will move simultaneously





❖ Align spectra (al) with each other by matching peak positions

- reference one spectrum to another
- the names of the axis have to be the same
- check the shift offset values in spectrum settings (st)

Overlay Spectra

❖ Overlay (ol) the contours of one spectrum onto another spectrum view

- spectrum type has to be the same (dimension & axis names)
- the contour levels of each spectrum are separately adjusted
- overlay as many spectra onto one view
- only the peaks and labels of the original view are shown




Peak Assignment

❖ Assign the selected peak using the assignment dialog (at)

- resonances known in the current project and closest in frequency to the peak position for the selected axis will be highlighted and can be selected
- For new resonances the residue group and atom name have to be entered manually
- rename, merge or swap resonances (rr)
- remove unused project resonances (dr)

❖ Assignment guessing

- assist peak assignment based on known resonances (select ppm ranges and the optional group and atoms)
- the values are also used by the assignment guessing pointer mode (sh-F2)




Lists of peaks and table of assignments

❖ Every spectrum in sparky has a peak list (lt) associated

- double clicking on a peak in this list will show and center this peak in the spectrum (or select the peak and type vc)
- output is customisable and exportable to other applications

❖ Assignment table (tb)


- multiple tables possible by defining different molecules or conditions in the spectrum settings (st)



Lists of resonances

❖ Resonance list (rl)

- convenient list of the combined resonance assignments
- inspect the standard deviation of each resonance to find assignment mistakes
- select the resonance to open the resonance assignments (ra) which lists all the peaks that have been assigned for this resonance - trace the deviating peak assignments and jump to the corresponding peak in the spectrum



More features in Sparky

❖ 3D, 4D, ..., spectra

- restricted peak picking
- strip plots

❖ Connection to structure

- view NOE crosspeak - Atom-Distance
- assignment assistance
- completeness check with spin graphs

❖ Relaxation Data

- fit the exponential decay from a series of 2D spectra

